陕西师范大学  
SHAANXI NORMAL UNIVERSITY

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# China Service Training Center --- Training Guide for Customer

Customer Enrollment Website (客户报名网址):

<http://www.thermo.com.cn/TrainingCenter>

On-line Training Video Entrance (在线视频自学网址):

<http://www.thermo.com.cn/TrainingCenter/VideoCourse-List.aspx>

On-line Training Registration (远程在线培训选课报名):

<http://www.thermo.com.cn/TrainingCenter/AirClassRoom.aspx>

Training Material Search and Download (资料下载专区):

<http://www.thermo.com.cn/TrainingCenter/DocumentSet.aspx>

联系我们: [training.china@thermofisher.com](mailto:training.china@thermofisher.com)/4006505118

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" ... the preparation obtained I name a *chromatogram* and the method proposed *chromatography* ... "

" ... it is obvious that the adsorption phenomena described are characteristic not only of chlorophyll pigments; it is evident that different coloured and colourless compounds follow the same regularities."

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
Properties Report Column

Categories:

Sequence  
Sample  
Audit Trail  
Preconditions  
Chromatogram  
Detection Parameters  
Peak Results  
Peak Calibration  
Peak Table  
Peak Purity and Identifi  
Quantification Method  
Information Table

Variables:

Signal Value at Peak End  
Baseline Value at Peak St  
Baseline Value at Peak Re  
Baseline Value at Peak En  
Detection Code at Peak St  
Detection Code at Peak En  
Type  
Modified  
Manually Assigned  
Resolution  
Asymmetry  
Theoretical Plates

OK

Cancel

Customize...

plain Variable...

Parameter...

Formula... peak asymmetry("ep")

Header... "Asym. "

Dimension... "(EP)"

Format... 0.00

Channel

Selected Channel

Fixed Channel (|IV\_VIS\_1

Parameter Input for 'Peak Asymmetry'

Compute 'Peak Asymmetry' according to the

EP-Formula (default): same as USP

AIA-Formula

Formula using Statistical Moments

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
Properties Report Column

Categories:

Sequence  
Sample  
Audit Trail  
Preconditions  
Chromatogram  
Detection Parameters  
Peak Results  
Peak Calibration  
Peak Table  
Peak Purity and Identifi  
Quantification Method  
Information Table

Variables:

Signal Value at Peak Rete  
Signal Value at Peak End  
Baseline Value at Peak St  
Baseline Value at Peak Re  
Baseline Value at Peak En  
Detection Code at Peak St  
Detection Code at Peak En  
Type  
Modified  
Manually Assigned  
Resolution  
Approximate

OK

Cancel

Customize...

Plain Variable

Formula... peak resolution ("ep")

Header... "Resol."

Dimension... "(EP)"

Format... 0.00

Channel

Selected Channel

Fixed Channel (UV\_VIS\_1)

Parameter Input for 'Peak Resolution'

Compute 'Peak Resolution' using

EP-Formula (default)

USP-Formula


Statistical Moments

The second peak in the resolution formula is

the next main peak in the chromatogram (

the previous main peak in the chromatogr

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